



Report for CoViDiag, Infected patient #1 (A1)

Sample: Infected patient #1
 Lot number:
 Well number: A1
 Configuration: CoViDiag 3.1 (1.0)

Image before analysis:



Analyzed image:



Replicate name	Result	Value
SARS-CoV-2		Intensity [AU]
Nucleocapside	Positive	77.290
Spike 1	Positive	65.437
Spike 2	Positive	83.957
Spike RBD	Positive	72.732
Spike NTD	Positive	78.757
		Intensity [AU]
Control +	Positive	73.523 VALID (min. 10.000 AU)
SampCtrl	Positive	65.527 VALID (min. 10.000 AU)

Test name	Result	Comment
SARS-CoV-2	Positive	


Test is valid.




Report for CoViDiag, Infected patient #2 (C1)

Sample: Infected patient #2
 Lot number:
 Well number: C1 †
 Configuration: CoViDiag 3.1 (1.0)

Image before analysis:



Analyzed image:



† grid changed by user

Replicate name	Result	Value
SARS-CoV-2		Intensity [AU]
Nucleocapside	Positive	70.969
Spike 1	Positive	30.812
Spike 2	Positive	54.844
Spike RBD	Positive	66.513
Spike NTD	Negative	0.973
		Intensity [AU]
Control +	Positive	60.093 VALID (min. 10.000 AU)
SampCtrl	Positive	61.003 VALID (min. 10.000 AU)

Test name	Result	Comment
SARS-CoV-2	Positive	

Test is valid.


Report for CoViDiag, Vaccinated patient (F2)

Sample: Vaccinated patient
 Lot number:
 Well number: F2
 Configuration: CoViDiag 3.1 (1.0)

Image before analysis:



Analyzed image:



Replicate name	Result	Value
SARS-CoV-2		Intensity [AU]
Nucleocapside	Negative	0.687
Spike 1	Positive	63.445
Spike 2	Positive	42.509
Spike RBD	Positive	76.773
Spike NTD	Negative	9.808
		Intensity [AU]
Control +	Positive	70.716 VALID (min. 10.000 AU)
SampCtrl	Positive	63.357 VALID (min. 10.000 AU)

Test name	Result	Comment
SARS-CoV-2	Positive	

Test is valid.




Report for CoViDiag, Non immunized patient (H2)

Sample: Non immunized patient
 Lot number:
 Well number: H2
 Configuration: CoViDiag 3.1 (1.0)

Image before analysis:



Analyzed image:



Replicate name	Result	Value
SARS-CoV-2		Intensity [AU]
Nucleocapside	Negative	0.576
Spike 1	Negative	0.439
Spike 2	Negative	0.218
Spike RBD	Negative	1.263
Spike NTD	Negative	0.776
		Intensity [AU]
Control +	Positive	72.902 VALID (min. 10.000 AU)
SampCtrl	Positive	78.107 VALID (min. 10.000 AU)

Test name	Result	Comment
SARS-CoV-2	Negative	

Test is valid.

